




THE NEURO LIGHT
- TOTAL BRAIN WELLNESS -

BASIC PROPOSITION



- Contrary to long held medical theory, the Human Brain is capable of dynamic change throughout our entire lifespan.
- This dramatic fact opens the doors for new avenues of increasing and sustaining Brain Wellness.

THE NEURO LIGHT: TOTAL BRAIN WELLNESS



- The Neuro Light is a multi-component device capable of positive effects in Brain Wellness related functions.
- The device is a combination of unique hardware and software designs.

THE NEURO LIGHT

Brain Gym

Neuroplasticity Training **Brain Food**

Transcranial Light Stimulation

Brain Skills


Function Specific Brain Training





- The theory of Neuroplasticity can be summed up as:
 - "If you FIRE it, you WIRE it".

- The brain responds to incoming stimulation as a source of information and guidance in the formation of novel patterns of adaptive behavior.
- Photonic (LIGHT) stimulation is a powerful and efficient type of brain stimulation.






- The <u>Neuro Light</u> recognizes that every individual has their own potential for Neuroplastic change.
- The potential relates to age, vitality and possible pathology.
- The term to describe this is **Neuroplastic Capacity**.

- The BRAIN GYM focuses on Neuroplastic Capacity.
- There are three sections in the BRAIN GYM:
 - 1) Neuroplastic Capacity ASSESSMENT
 - 2) Neuroplastic Capacity TRAINING
 - 3) Network Reinforcements





- There are three ways that Neuroplastic Capacity can increase:
 - 1) Synaptic Plasticity new pathways using pre-existing neurons
 - 2) Neurogenesis new neurons to serve new pathways
 - 3) Functional Compensatory Plasticity new pathways are formed in areas of the brain that do not normally service these functions (especially important in aging, brain trauma and diseased brains)



- The Neuro Light has three (3) levels of Neuroplasticity Training:
 - Foundation
 - Intermediate
 - Advanced

An individual progresses through the levels of training in stages.



- Each level type of Neuroplastic Capacity <u>Training</u> has elements designed to increase <u>Neuroplastic</u> Capacity in an integrated dynamic.
- The <u>analogy</u> of physical exercise is helpful.

- Physical exercise analogy as related to the Neuroplasticity Training in the Neuro Light:
- Coordination > Signal <u>Frequency</u>
- Flexibility > Signal <u>Pattern</u>
- Endurance > Signal <u>Change</u>
- Strength > Signal Intensity





- These three (3) levels of Neuroplastic Capacity (Foundation, Intermediate, Advanced) are basic concepts to help in understanding the application of the Neuro Light.
- No one knows exactly how to measure the Neuroplastic Capacity of a person.
- The best you can do is make an Assessment which approximates the Wellness of the brain.



- The results of a Neuroplastic Capacity <u>Assessment</u> allow you to:
 - 1) focus on good choices of specific Neuroplastic Capacity Trainings.
 - 2) have guidance in what Network Reinforcements and/or Brain Skills Trainings to choose.
 - 3) All Brain Skills Trainings are coded to match the levels of Neuroplastic Capacity.



- There are 3 ways to do a Neuroplastic Capacity Assessment:
 - 1) Subjective Assessment
 - 2) Q & A Interview Assessment
 - 3) Neuro Light Experience Assessment



- Keep in mind that the Neuroplastic Capacity of a person is complex.
- No assessment is 100% perfect.
- Results of any assessment should be used as broad guide subject to reassessment based on real experience.



- #1 Subjective Assessment:
 - These are your basic impressions of the person.
 - Use the analogy of physical exercise as a basis for your assessment.



- #1 Subjective Assessment:
 - Because of the Mind-Body relationship, more than just an analogy,
 you can use their physical body to help inform your Subjective Assessment.
 - It is well known that physical exercise and activities have a strong correlation to levels of Neuroplastic Capacity.



- #1 <u>Subjective</u> <u>Assessment</u> (<u>examples</u>)
 - Coordination are they mentally smooth and graceful?
 - Flexibility do they seem open minded and fluid in thought?
 - Endurance can they maintain a higher level of communication over time?
 - Strength are they able to project intensity with balanced behavior?



- #1 Subjective Assessment
- Rate each category at a level of Neuroplastic Capacity (example):
 - Coordination Foundation
 - Flexibility Intermediate
 - Endurance Intermediate
 - Strength Foundation



- #1 Subjective Assessment
- OK. Now based on your <u>Subjective Assessment</u> you know:
 - 1) In the BRAIN GYM, what Neuroplastic Trainings will be helpful.
 - 2) in the BRAIN SKILLS, what Brain Training sessions will likely to be appropriate.




- #2 Q & A Interview Assessment:
 - This style of assessment can stand on its own or supplement your Subjective Assessment.
 - The Q & A can be helpful in clarifying areas where you may not feel so confident in your Subjective Assessment.



- #2 Q & A Interview Assessment (examples)
 - Coordination how well can you balance and entertain multiple thoughts?
 - Flexibility can you easily move in and out of different emotions and ideas?
 - Endurance is it easy for you to keep your focus and concentration?
 - Strength do you feel you have the mental power and memory you need?

- #2 Q & A Interview Assessment
- Similar to the Subjective Assessment...
 <u>rate</u> each category at a level of Neuroplastic Capacity (<u>example</u>):
 - Coordination Advanced
 - Flexibility Intermediate
 - Endurance Intermediate
 - Strength Foundation







- #3 NEURO LIGHT Experience Assessment:
 - In the BRAIN GYM section of the Neuro Light, there are specific Neuroplastic Capacity Assessment processes.
 - These Light Experiences are another way of helping you to make an assessment of a person.



- #3 NEURO LIGHT Experience Assessment:
 - Each of the 4 aspects of basic Neuroplastic Capacity have their own Light Experience assessment.
 - Coordination, Flexibility, Endurance, Strength

- #3 NEURO LIGHT Experience Assessment:
 - The <u>structure</u> of each <u>Light Experience</u> is the <u>same</u>:
 - A) it begins with a short period of darkness
 - B) then a Light Experience of about 30 40 seconds
 - C) another period of darkness for about 20 seconds
 - D) a second Light Experience of about 30 40 seconds
 - E) another period of darkness for about 20 seconds
 - F) a third Light Experience of about 30 40 seconds
 - G) a final period of darkness





- #3 NEURO LIGHT Experience Assessment:
- Before the Light Experience, the person is told that there will be three short Light Experiences separated by short periods of darkness.
- They are told that there will be one of the three that is "most comfortable".
- After the three Light Experiences, they will choose one as "most comfortable".




- #3 NEURO LIGHT Experience Assessment:
 - The sequence of the 3 Light Experiences is always the same.

(Intermediate, Foundation, Advanced)

 Do <u>not</u> tell the person in advance what the levels are in the <u>Light Experiences</u>



- #3 NEURO LIGHT Experience Assessment:
 - Based on the choices made by the person of the "most comfortable", you can assign a level (Foundation, Intermediate, Advanced) of Neuroplastic Capacity in that particular area (Coordination, Flexibility, Endurance, Strength).




- Having performed one or more Neuroplastic Capacity <u>Assessments</u>, you now know where to focus the following steps:
 - 1) what Neuroplastic Capacity Trainings to use
 - 2) what Brain Skills Trainings to use



- Keep in mind that the results of the Assessments are general guidelines and NOT strict rules.
- The combination of Subjective, Q & A Interview and NEURO Light Experience Assessments will give a good general impression.
- Be creative and attempt to achieve the most while doing the least.



- Following your Neuroplastic Capacity Assessment(s), a typical session could be:
- 1) one (or two) 11 minute Neuroplastic Capacity Trainings
 Choose the one (or two) that appear to benefit from fortification
- 2) one 11 minute Network Reinforcement or Brain Skills Trainings
 - Choose the one (or two) that match the Assessment level of Capacity and that correspond to their general goals



- It is better to approach the NEURO LIGHT processes as an experience of EXERCISE & LEARNING.
- They are **NOT** medical treatments or therapies.

- After finishing the various Light Trainings, do a Summary Discussion with the person.
- Ask them about their feelings and experiences and LISTEN to what they say.
- Do NOT project your ideas and concepts onto them.

- The Mirror approach works very well:
- 1) Listen to what they say
- 2) Pick out the "hot word" in their statement
- 3) Feedback that "hot word" in the form of a question or just repeat it
- 4) They WILL respond further based on the "hot word"
- 5) Repeat this same cycle 4 or 5 times with new "hot words"
- 6) An appreciation or understanding will typically present itself

Strategy without action is worthless.

- Final Step convert the "understanding" into an Action Strategy
- 1) They have a basic reason or goal in mind which is why they are having Neuro Light sessions
- 2) Ask them if the "understanding" that has emerged from the session has any possible correlation with their reason or goal
- 3) Ask them to commit to some (small?) Action in the following days that may help in integrating the session experience into their actual daily life.

- The Neuro Light is capable of generating select types of Photonic signals related to Frequency, Pattern, Change, Intensity.
- These Neuroplastic "exercises" result in increased capacity for higher levels of <u>FIRE it – WIRE it</u> responses in the brain.



BRAIN SKILLS FUNCTION SPECIFIC BRAIN TRAINING


Remember that for easy and accurate selection, the 80 Function Specific Brain Training programs in the Brain Skills are subdivided to correspond to the same levels found in the Brain Gym (Neuroplasticity Trainings)

- Foundation
- Intermediate
- Advanced



BRAIN SKILLS FUNCTION SPECIFIC BRAIN TRAINING

- For further ease of selection, the Brain Training sessions in Brain Skills are also in general ZONES:
 - Peace
 - Cognition
 - Mood
 - Dependency
 - Exploration



THE NEURO LIGHT - FULL INTEGRATION - TOTAL BRAIN WELLNESS

Brain Gym

Neuroplasticity Training **Brain Food**

Transcranial Light Stimulation

Brain Skills

Function Specific Brain Training

